

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HAB01-01	Crude antisera	2 rabbits	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	60-80 ml crude antisera Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HAB01-02	Crude antisera	2 rabbits	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; ELISA, guaranteed titer >1:10 000	100 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HAB01-03	Crude antisera	2 rabbits	Package 3	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 6 immunisations, 90 days ; ELISA, guaranteed titer >1:10 000	200 ml crude antisera Pre-immune sera Report with ELISA data 10 mg free peptide, purity >90%, LC & MS analysis	4.5 month
HAB02-01	Crude antisera	2 guinea pigs	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	5-8 ml crude antisera Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HAB02-02	Crude antisera	2 guinea pigs	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; ELISA, guaranteed titer >1:10 000	10-15 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HAB03-01	Crude antisera	2 rats	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	2-3 ml crude antisera Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HAB03-02	Crude antisera	2 rats	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; ELISA, guaranteed titer >1:10 000	5-8 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HAB04-01	Crude antisera	2 mice	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	0.5-1 ml crude antisera Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HAB04-02	Crude antisera	4 mice	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	1-2 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HAB05-01	Crude antisera	2 Chickens	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; ELISA, guaranteed titer >1:10 000	5-10 ml crude antisera Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3 month
HAB06-01	Crude antisera	1 Goat	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 6-8 immunisations, 150 days ; ELISA, guaranteed titer >1:10 000	500 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	5 month
HAB06-02	Crude antisera	1 Goat	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 10-12 immunisations, 245 days ; ELISA, guaranteed titer >1:10 000	1000 ml crude antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	9 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP1AB01-01	Total IgG purified	2 rabbits	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	15-25 mg total IgG Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP1AB01-02	Total IgG purified	2 rabbits	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	100-200 mg total IgG Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HP1AB01-03	Total IgG purified	2 rabbits	Package 3	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 6 immunisations, 90 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	300-400 mg total IgG Pre-immune sera Report with ELISA data 10 mg free peptide, purity >90%, LC & MS analysis	4.5 month
HP1AB02-01	Total IgG purified	2 guinea pigs	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	5-10 mg total IgG Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP1AB02-01	Total IgG purified	2 guinea pigs	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	30-50 mg total IgG Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HP1AB03-01	Total IgG purified	2 rats	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	5-10 mg total IgG Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3 month
HP1AB03-02	Total IgG purified	2 rats	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	20-40 mg total IgG Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP1AB04-01	Total IgG purified	2 mice	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	4-5 mg total IgG Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3 month
HP1AB04-02	Total IgG purified	4 mice	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Chromatography Protein A ; Quantitation ; ELISA, guaranteed titer >1:20 000	8-10 mg total IgG Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP1AB05-01	Total IgY purified	2 Hens	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4 immunisations, 70 days ; Egg collection, yolk isolation and total IgY purification ; Quantitation ; ELISA, guaranteed titer >1:20 000	50-100 mg total IgY Pre-immune purified total IgY Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP1AB05-02	Total IgY purified	2 Hens	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Egg collection, yolk isolation and total IgY purification ; Quantitation ; ELISA, guaranteed titer >1:20 000	200-300 mg total IgY Pre-immune purified total IgY Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HP1AB06-01	Total IgG purified	1 Goat	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 6-8 immunisations, 150 days ; Chromatography Protein G ; Quantitation ; ELISA, guaranteed titer >1:20 000	Purification of 500 ml antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	5 month
HP1AB06-02	Total IgG purified	1 Goat	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 10-12 immunisations, 245 days ; Chromatography Protein G ; Quantitation ; ELISA, guaranteed titer >1:20 000	Purification of 1 000 ml antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	9 month
Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP2AB01-01	Affinity purified	2 rabbits	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4 immunisations, 70 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	2-5 mg specific antibodies Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP2AB01-02	Affinity purified	2 rabbits	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4-6 immunisations, 70-90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	6-10 mg specific antibodies Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month
HP2AB01-03	Affinity purified	2 rabbits	Package 3	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;6 immunisations, 90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	11-20 mg specific antibodies Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4.5 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP2AB02-01	Affinity purified	2 guinea pigs	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4 immunisations, 70 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	2-3 mg specific antibodies Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP2AB02-01	Affinity purified	2 guinea pigs	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4-6 immunisations, 70-90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	4-10 mg specific antibodies Pre-immune sera Report with ELISA data 2-3 mg free peptide, purity >90%, LC & MS analysis	4 month
HP2AB03-01	Affinity purified	2 rats	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	0.2 - 1 mg specific antibodies Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP2AB04-01	Affinity purified	2 mice	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	50-100 µg specific antibodies Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP2AB04-02	Affinity purified	4 mice	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 4-6 immunisations, 70-90 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	100-200 µg specific antibodies Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	3.5 month
HP2AB05-01	Affinity purified	2 Hens	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4 immunisations, 70 days ; Egg collection, yolk isolation and total IgY purification ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	5-10 mg specific antibodies Pre-immune purified total IgY Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4 month

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP2AB05-02	Affinity purified	2 Hens	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;4-6 immunisations, 70-90 days ; Egg collection, yolk isolation and total IgY purification ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	11-20 mg specific antibodies Pre-immune purified total IgY Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	4.5 month
HP2AB06-01	Affinity purified	1 Goat	Package 1	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ; 6-8 immunisations, 150 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	Purification of 500 ml antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	5.5 month
HP2AB06-02	Affinity purified	1 Goat	Package 2	10-15 aa peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your peptide ;10-12 immunisations, 245 days ; Preparation of affinity column by covalent linkage of your peptide ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	Purification of 1 000 ml antisera Pre-immune sera Report with ELISA data 5 mg free peptide, purity >90%, LC & MS analysis	9.5 month
Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime

Antikörperherstellung - Haptene: Synthetische Peptide

Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HAB07-01	Peptide synthesis for immunisation		1 aa		Additional peptide synthesis cost per aa for peptides >15 aa	n/a	3-5 days
HAB07-02	Antigen conjugation (BSA or OVA)				Conjugation to BSA or OVA of your peptide and purification of the conjugate for immunisation	n/a	3-5 days
HAB07-03	Peptide purity > 95%				Purification > 95% of your peptide	n/a	3-5 days
HAB07-04	Multipptide Antibodies	All	1 additional peptide	10-15 aa additional peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your additional peptide ; ELISA against additional peptide	2-5 mg of additional peptide	
Ref	Type	Host	Size	Hapten	Service features	Delivery	Leadtime
HP2AB08-01	Phospho-specific antibodies	2 rabbits	Package 1	10-15 aa phospho-peptide with terminal cys for conjugation	Protein analysis and peptide design ; Synthesis, purification & KLH conjugation of your phospho-peptide and non-phosphopeptide; 4 immunisations, 70 days ; Preparation of two affinity column by covalent linkage of your peptides ; Affinity chromatography of your antibodies ; Quantitation ; ELISA, guaranteed titer >1:20 000	2-5 mg phospho-specific antibodies 3-5 mg non-phospho-specific antibodies Pre-immune sera Report with ELISA data and cross-reaction analysis 2-3 mg free peptides, purity >90%, LC & MS analysis	4 month