

Cat.No	Name	Sensitivity	Detection Range	Reaction Type
T-SEKH-0001	Human IL-1 α ELISA Kit	4 pg/ml	7.81-500 pg/ml	Sandwich-ELISA
T-SEKH-0002	Human IL-1 β ELISA Kit	7.8 pg/ml	15.62-1000 pg/ml	Sandwich-ELISA
T-SEKH-0002(HS)	Human IL-1 β (High Sensitivity) ELISA Kit	0.15 pg/ml	0.312-20 pg/ml	Sandwich-ELISA
T-SEKH-0003	Human IL-1F9 ELISA Kit	18.75 pg/ml	37.5-2400 pg/ml	Sandwich-ELISA
T-SEKH-0004	Human IL-1RA ELISA Kit	39 pg/ml	78.125-5000 pg/ml	Sandwich-ELISA
T-SEKH-0005	Human IL-1 RI ELISA Kit	62.5 pg/ml	125-8000 pg/ml	Sandwich-ELISA
T-SEKH-0006	Human IL-1 RII ELISA Kit	15 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA
T-SEKH-0007	Human IL-1R4 ELISA Kit	15 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA
T-SEKH-0008	Human IL-2 ELISA Kit	3.9 pg/ml	7.81-500 pg/ml	Sandwich-ELISA
T-SEKH-0008(HS)	Human IL-2(High Sensitivity) ELISA Kit	0.3 pg/ml	0.625-40pg/ml	Sandwich-ELISA
T-SEKH-0009	Human IL-2R α ELISA Kit	31.25 pg/ml	62.5-4000 pg/ml	Sandwich-ELISA
T-SEKH-0010	Human IL-3 ELISA Kit	7 pg/ml	15.625-1000 pg/ml	Sandwich-ELISA
T-SEKH-0011	Human IL-4 ELISA Kit	2 pg/ml	3.12-200 pg/ml	Sandwich-ELISA
T-SEKH-0011(HS)	Human IL-4(High Sensitivity) ELISA Kit	0.15 pg/ml	0.625-40pg/ml	Sandwich-ELISA
T-SEKH-0012	Human IL-5 ELISA Kit	2 pg/ml	4.68-300 pg/ml	Sandwich-ELISA
T-SEKH-0013	Human IL-6 ELISA Kit	2 pg/ml	3.9-250 pg/ml	Sandwich-ELISA
T-SEKH-0013(HS)	Human IL-6(High Sensitivity) ELISA Kit	0.3 pg/ml	0.625-40 pg/ml	Sandwich-ELISA
T-SEKH-0014	Human IL-6R α ELISA Kit	4 pg/ml	7.81-500 pg/ml	Sandwich-ELISA
T-SEKH-0015	Human IL-7 ELISA Kit	3.9 pg/ml	7.81-500 pg/ml	Sandwich-ELISA
T-SEKH-0016	Human IL-8 ELISA Kit	15 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA
T-SEKH-0016(OS)	Human IL-8(One-Step) ELISA Kit	15.625 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA