

Cat.No	Name	Sensitivity	Detection Range	Reaction Type	Assay Time
T-SEKM-0031(OS)	Mouse IFN- γ (One-Step) ELISA Kit	15.625pg/ml	31.25-2000pg/ml	Sandwich-ELISA	60min
T-SEKM-0032	Mouse IFN- β ELISA Kit	12.5 pg/ml	25-1600 pg/ml	Sandwich-ELISA	3h
T-SEKM-0032(OS)	Mouse IFN- β (One-Step) ELISA Kit	12.5pg/ml	25-1600pg/ml	Sandwich-ELISA	60min
T-SEKM-0033	Mouse IFN-zeta ELISA Kit	7.81 pg/ml	15.63-1000 pg/ml	Sandwich-ELISA	3h
T-SEKM-0034	Mouse TNF- α ELISA Kit	15.63 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	3h
T-SEKM-0034(HS)	Mouse TNF- α (HS) ELISA Kit	0.78 pg/ml	1.56-100 pg/ml	Sandwich-ELISA	3h
T-SEKM-0034(OS)	Mouse TNF- α (One-Step) ELISA Kit	15.625 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	60min
T-SEKM-0035	Mouse TGF- β 1 ELISA Kit	3.9 pg/ml	7.81-500 pg/ml	Sandwich-ELISA	3h
T-SEKM-0036	Mouse TGF- β 2 ELISA Kit	1.95 pg/ml	3.91-250 pg/ml	Sandwich-ELISA	3h
T-SEKM-0037	Mouse VCAM-1 ELISA Kit	15.63 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	3h
T-SEKM-0038	Mouse EGF ELISA Kit	1.95 pg/ml	3.9-250 pg/ml	Sandwich-ELISA	3h
T-SEKM-0039	Mouse VEGF ELISA Kit	7.81 pg/ml	15.63-1000 pg/ml	Sandwich-ELISA	3h
T-SEKM-0040	Mouse G-CSF ELISA Kit	7.81 pg/ml	15.63-1000 pg/ml	Sandwich-ELISA	3h
T-SEKR-0055	Rat ANG1 ELISA Kit	12 pg/ml	23.43-1500 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0056	Rat ANG2 ELISA Kit	40 pg/ml	81.25-5200 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0057	Rat IL-12 ELISA Kit	31.25 pg/ml	62.5-4000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0058	Rat BNP ELISA Kit	0.08 ng/ml	0.15-10 ng/ml	Sandwich-ELISA	3.5h
T-SEKR-0059	Rat CK-MB1 ELISA Kit	250 pg/ml	1.25-80 ng/ml	Sandwich-ELISA	3.5h
T-SEKR-0060	Rat AHSR ELISA Kit	9.3 ng/ml	18.75-1200 ng/ml	Sandwich-ELISA	3.5h
T-SEKR-0061	Rat MIF ELISA Kit	15.625 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0062	Rat HGF ELISA Kit	39 pg/ml	78.12-5000 pg/ml	Sandwich-ELISA	1.3h
T-SEKR-0063	Rat Adiponectin ELISA Kit	21.75 pg/ml	43.75-2800 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0064	Rat Prolactin ELISA Kit	46.8 pg/ml	93.75-6000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0065	Rat MMP-1 ELISA Kit	93.75 pg/ml	187.5-12000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0066	Rat MMP-2 ELISA Kit	80 pg/ml	0.31-20 ng/ml	Sandwich-ELISA	3.5h

T-SEKR-0067	Rat MMP-3 ELISA Kit	12.5 pg/ml	25-1600 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0068	Rat MMP-7 ELISA Kit	25 pg/ml	50-3200 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0069	Rat MMP-8 ELISA Kit	46 pg/ml	93.75-6000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0070	Rat MMP-14 ELISA Kit	125 pg/ml	250-16000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0072	Rat IL-23 ELISA Kit	31.25 pg/ml	62.5-4000 pg/ml	Sandwich-ELISA	3.5h
T-SEKR-0073	Rat MPO ELISA Kit	1.6 ng/ml	3.125-200 ng/ml	Sandwich-ELISA	3.5h
T-SEKRT-0004	Rabbit IL-8 ELISA Kit	9 pg/ml	18.75-1200 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0006	Rabbit IL-5 ELISA Kit	15 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0009	Rabbit VEGF ELISA Kit	37 pg/ml	75-4800 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0019	Rabbit IL-1 β ELISA Kit	15 pg/ml	31.25-2000 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0020	Rabbit IL-6 ELISA Kit	17 pg/ml	34.37-2200 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0023	Rabbit PLA2G7 ELISA Kit	200 pg/ml	468.75-30000 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0027	Rabbit HGF ELISA Kit	34 pg/ml	68.75-4400 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0030	Rabbit PDGF-AB ELISA Kit	9 pg/ml	18.75-1200 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0032	Rabbit PIGF ELISA Kit	3 pg/ml	6.25-400 pg/ml	Sandwich-ELISA	2.33h
T-SEKRT-0033	Rabbit IL-2 ELISA Kit	23 pg/ml	46.87-3000 pg/ml	Sandwich-ELISA	2.33h

Send your request at:

Email: info@sansunbiologics.com